




Fetal Scalp Blood Sampling Procedure Pack

Reduce risk with the BridgeMaster™ Fetal Blood Sampling Wand.

- Single handed sampling wand with a maximum incision depth of 1.7mm.
- Wand is graduated for assessment of fetal descent.
- No petroleum jelly required.
- Spray-coated heparinised capillary tubes and mixing wires help to prevent clotting.
- Pre-mounted retractable blade prevents fingerstick injury and reduces preparation time.
- Capillary tube adapter provided is compatible with the majority of blood gas analysers.
- Five capillary tubes and caps included for multiple analyses in the same pregnancy.
- A sterile drape is provided.



NICE Clinical Guideline 132 Recommendation 1.3.2.5

The ease of taking a fetal blood sample is crucial in providing reliable pH results, as they are to be based on the 'gold standard' cut off values in Table 1.

рН	Description	Action	
>7.2	Normal fetus	Blood sampling should be repeated if the fetal heart rate abnormality persists	
7.2 - 7.25	Pre-acidotic	Repeat blood sample within 30 minutes or consider delivery if rapid fall since last sample	
<7.2	Acidotic	Acidotic Delivery immediately	

Table 1: Cut off values of fetal blood pH and the actions to take.



Midwife





Fetal scalp blood sampling

ELA01	Electrically Illuminated Amnioscope (28mm)	5 Single Use Units per Case
ELA02	Omniscope (18mm)	5 Single Use Units per Case
FBS-AC	Fetal Scalp Blood Sampling Procedure Pack	Includes: 1 x BridgeMaster™ FBS Fetal Blood Sampler/6 x Long Cotton Tipped Buds/5 x Heparinised Capillary Tubes/ 10 x Capillary Tube End Caps/ Mixing Wires/ 4 x Universal Blood Gas Analyser Adapters with clot buster/ 1x Sterile Surgical Drape

Cruinn is proud of its commitment to provide excellent customer service and support; with the following features:

- Online Training Programs
- 24/7 UK Technical and Clinical Support
- Guaranteed free of charge replacement devices provided in the event of technical issues.

Very easy to use, and much more comfortable to the patient than the conventional.

R.Singh, Specialist Registrar

