



Unique monoclonal
and polyclonal antibody
approach unites to create the
most sensitive ROM test.

Shh... we're hunting markers.

ROM
plus⁺

Rupture of Membranes Test



Accurate diagnosis when it matters most

- ROM Plus ES demonstrated 100% sensitivity and 100% specificity between 24-34 weeks gestation¹

"...development of a point-of-care test that is both accurate and easy to use by a wide array of clinicians (nurses, midwives, physicians, etc.) would be a true asset."

—Mariona F. Cabero L. J Matern Fetal Neonatal Med. 2012.

ROM⁺
plus[®]

Rupture of Membranes Test



Jan Bye a.s

MEDISINSK TEKNIKK

www.janbye.no post@janbye.no

Ordering Information


Box

POC (Point of Care) Test Kits ROM Plus® ES Complete Test Kit Test kit includes a lateral flow test strip, swab, vial with buffer solution and convenient shipper cap	ROM-1025	25 ea
ROM Plus® ES Strip Quality Control 1 positive control and 1 negative control. Each box contains 1 activation sleeve	ROM-QCS	5 ea



1. Thomasino T, Levi C, Draper M, Neubert AG. Diagnosing rupture of membranes using combination monoclonal/polyclonal immunologic protein detection. *J Reprod Med.* 2013;58(5-6):187-94.
2. Mariona FG, Cabero L. Are we ready for a new look at the diagnosis of premature rupture of membranes? *J Matern Fetal Neonatal Med.* 2012;25(4):403-407.
3. Senanayake HM. Actim™ PROM, AmniSure®, and ROM+plus®: Rupture of membrane kits tested on amniotic fluid from women at C-section: a comparative study. *Sri Lanka Journal of Obstetrics and Gynaecology.* 2013; 116-121. Available at: <http://www.slcoog.lk/sljog/>

Your Local Distributor:




Jan Bye a.s

MEDISINSK TEKNIKK

www.janbye.no post@janbye.no


Test Procedure

Collect sample with vaginal swab




1

Mix swab in buffer vial




2

Add test strip to vial and set timer



3

Read test result



4



**CLINICAL
INNOVATIONS
EUROPE**

6-9 The Square
Stockley Park
Heathrow

United Kingdom
UB11 1FW

Phone: +44 (0)333-370-4408
Fax: +44 (0)845-299-1517
www.clinicalinnovations.com

Not made with natural rubber,
latex, PVC or DEHP.

Made in the USA

Patents Pending

© Copyright 8-2014



Rupture of Membranes Test